

NEPTUNE

Liquid filling table tops

Cleaning

- Exceptional cleaning process of our system, resulting in grater operational efficiency.

Dosage range

- Volumes ranging from 4 to 750 ml, depending on specific requirements. Possibility of using up to two dosers working in parallel.

Precision

- Dosing accuracy well below 1% for the above-mentioned dosages.

Versatility

- Used to fill liquids in multiple packages, it can precisely handle a wide variety of fluids regardless of their density and viscosity.

Competitive pricing

- Available at a lower price than conventional dosing systems.

Liquid dosing technologies



**Volumetric
Dosing for Water-
Based Liquids**



**Volumetric
Dosing for
Dense Fluids**



**Volumetric
Dosing for
Large Volumes**



**Volumetric
Dosing for
Dense Fluids**

Flow Controlling



**Compression
valve**



**Diaphragm
valve**

- Valve body material may vary depending on the target sector: Nylon PA12, Anticorodal-100, stainless steel 316L. All these valves are pneumatically operated.