

ZEPHYR

Powder dispenser



Dosage

- Doses from 50 mg to 2000 mg

Productivity

- 150 vpm

Disk diameter


- 200 mm

Disk format

- From 5 to 10 formats

Exceptional accuracy

- Up to 5% dosing error rate




Easy cleaning

- Easy washable product contact materials, ensuring high operation standards

Intuitive design

- Air/vacuum passage channels of larger dimensions than competitors, improved core formation fluid dynamics and ejection of the same with less turbulent motions



Total control

- Manual assembly using a PEEK ferrule
- Manual adjustment of the metered dose volume using a vernier
- In alternative to manual adjustment, possibility of the automatic adjustment of the volume compliance with the self-regulation for 100% in line check on dosing system

Services



**12 months
warranty**

Products



Granules



Powders

SETIN Equipment

Via Palazzetti, 5B, Int. 6, 40068,
San Lazzaro di Savena, Bologna, Italia

+39 0510019507

info@setinitalia.com



setinitalia.com